

Rigi10™ malleable penile prosthesis

Product Overview

For the Physicians

♂ Easy to implant

High bending angle may allow implantation with a smaller incision.

♂ Easy to size

Available in five diameters to better meet the patient's anatomical needs.

♂ Easy to assemble

SecureFit™ Extenders connect easily to add length. All diameters come with two each of the extender sizes (0.5 cm and 1.0 cm).

For the Patients

♂ Rigid shaft and distal end for a strong erection.

♂ High bending angle and memory for easy and comfortable concealment.

♂ Durable. Rigi10™ continues to offer satisfactory performance results even after long cycles of simulated tests.¹



Flexible Rod Technology™

Flexible Rod Technology™ provides easy bending and avoids undesirable turning on the axis of the implant when an erection is desired by the patient.

Flexible Rod Technology™ enables the rod to reach a bending angle of more than 135 degrees while simulating a more natural erection feeling.


SecureFit™ Extenders

Each Rigi10™ malleable penile prosthesis comes with a set of SecureFit™ Extenders, available in 0.5 cm and the 1.0 cm sizes.

SecureFit™ Extenders are easily attached to the Flexible Rod. The firm connection provides the physician with confidence in a secure connection between the extenders and the Flexible Rod.



Dimensional Specifications & Ordering Information



Prosthesis Product	Prosthesis Diameters	Prosthesis Lengths	Extenders
RG1009H	9 mm	23 cm	0.5 cm and 1 cm
RG1010H	10 mm	23 cm	
RG1011H	11 mm	25 cm	
RG1012H	12 mm	25 cm	
RG1013H	13 mm	25 cm	
RG1014H	14 mm	25 cm	

Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.

Prior to use, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.

References: 1 Data on file RIGI10-DOF-01-01.2019

Rigicon, INC.
2805 Veterans Memorial Hwy STE 5,
Ronkonkoma, NY 11779 United States of America
+1 (888) 202-9790 - www.rigicon.com

Rigicon®
Innovative Urological Solutions