

Rigi10™

malleable penile prosthesis

Product Overview

For the Physicians

♂ Easy to implant

High bending angle may allow implantation with a smaller incision. The increased lubricity with **hydrophilic coating** may promote easier surgical placement.

♂ Easy to size

Available in six diameters to better meet the patient's anatomical needs.

♂ Easy to assemble

SecureFit™ Extenders connect easily to add length. All diameters come with two each of the extender sizes (0.5 cm and 1.0 cm).

♂ Easy to prevent from infection

The hydrophilic coating offers physicians the freedom to choose their preferred aqueous solution.

For the Patients

♂ Rigid shaft and distal end for a strong erection.

♂ High bending angle and memory for easy and comfortable concealment.

♂ Durable. Rigi10™ malleable penile prosthesis continues to offer satisfactory performance results even after long cycles of simulated tests.¹



Flexible Rod Technology™

Flexible Rod Technology™ provides easy bending and avoids undesirable turning on the axis of the implant when an erection is desired by the patient.

Flexible Rod Technology™ enables the rod to reach a bending angle of more than 135 degrees while simulating a more natural erection feeling.

SecureFit™ Extenders

Each Rigi10™ hydrophilic malleable penile prosthesis comes with a set of SecureFit™ Extenders, available in 0.5 cm and the 1.0 cm sizes.

SecureFit™ Extenders are easily attached to the Flexible Rod. The firm connection provides the physician with confidence in a secure connection between the extenders and the Flexible Rod.




Hydrophilic Coating

Rigi10™ malleable penile prosthesis incorporates a hydrophilic coating on all external component surfaces. The hydrophilic coating offers you the freedom to choose the aqueous solution. Hydrophilic coating facilitates the rapid and strong absorption of the solution on the device and may promote easier device implantation.

Offers 6 different diameters

Diverse diameter options with Rigi10™ ensure best fit for the anatomical needs of each patient. Rigi10™ 9 mm and 10 mm rods are designed with conical distal tips to facilitate implantation in patients with narrower or compromised corpora cavernosa.



Ordering Information

Prosthesis Product	Prosthesis Diameters	Prosthesis Lengths	Extenders
RG1009H	9 mm	23 cm	0.5 cm and 1 cm
RG1010H	10 mm	23 cm	
RG1011H	11 mm	25 cm	
RG1012H	12 mm	25 cm	
RG1013H	13 mm	25 cm	
RG1014H	14 mm	25 cm	

Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.

Prior to use, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.

References: 1 Data on file RIGI10H-DOF-01-01.2019

Rigicon, INC.
2805 Veterans Highway Suite 5, Ronkonkoma,
NY 11779 United States of America
+1 (888) 202-9790 - www.rigicon.com

Rigicon®
Innovative Urological Solutions