

# Testi10™


## Testicular Prosthesis

Designed with a soft texture to offer patients missing one or both testicles a more natural feel

- » Available in five sizes to meet the patient's anatomical and sizing needs
- » Available pre-filled for simplified surgical procedure



## ■ Specifications and Ordering Information



Testi10™ Saline Pre-filled Testicular Prosthesis				PRE-FILLED
Product Code	Size	Volume	A (cm)	B (cm)
TESTIS-XS	Extra small	6 ml	2.7	2.1
TESTIS-S	Small	12 ml	3.2	2.6
TESTIS-M	Medium	18 ml	4.1	2.9
TESTIS-L	Large	22 ml	4.7	3.2
TESTIS-XL	Extra Large	26 ml	5.1	3.2

Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.

Refer to Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.

Rigicon, INC.  
2805 Veterans Highway Suite 5, Ronkonkoma,  
NY 11779 United States of America  
+1 (888) 202-9790 - [www.rigicon.com](http://www.rigicon.com)

**Rigicon**®  
Innovative Urological Solutions

TESTI10-20-11-2022