

Testi10™


Firm Testicular Prosthesis

Testi10™ is designed using high quality silicone, to help patients experience a complete and natural feeling.

- » Available in five sizes to meet the patient's anatomical and sizing needs



■ Specifications and Ordering Information



Testi10™ Firm Testicular Prosthesis				
Product Code	Size	Volume	A (cm)	B (cm)
TESTIF-XS	Extra small	6 ml	2.7	2.1
TESTIF-S	Small	12 ml	3.2	2.6
TESTIF-M	Medium	18 ml	4.1	2.9
TESTIF-L	Large	22 ml	4.7	3.2
TESTIF-XL	Extra Large	26 ml	5.1	3.2

Caution: Federal Law (U.S.) restricts this device to sale by or on the order of a physician.

Refer to Instructions for Use for a complete listing of indications, contraindications, warnings, precautions and potential adverse events.

Rigicon, INC.
2805 Veterans Highway Suite 5, Ronkonkoma,
NY 11779 United States of America
+1 (888) 202-9790 - www.rigicon.com

Rigicon®
Innovative Urological Solutions

TESTI0F-01-11-2022